

# BeVision M1

Move Beyond Vision

PARTICLE SIZE

PARTICLE SHAPE



## BeVision M1: Move Beyond Vision

The BeVision M1 provides an accurate analysis of particle size and shape in the range of 1 - 10,000  $\mu\text{m}$ . Besides, the BeVision M1 can be a vital part of the surface cleanliness test.

With the precise auto-moving stage, the BeVision M1 can do analysis more efficiently and go beyond the size limit.

The BeVision software helps you evaluate particle size and shape from 24 different aspects, and further organizes the data into an all-around validation of particles.

### Features and Benefits

- Measurement range: **1 - 10,000  $\mu\text{m}$**
- Results in compliance with **ISO 9276 - 6**
- Highly **reproducible** measurements
- A **high-speed** CCD camera
- **Powerful software** eases your work
- Automatic sample stage with **high position accuracy**
- **24** different particle size and shape parameters
- A key part of **surface cleanliness measurements**

With high magnification up to

**800** times\*

*\*Includes digital magnification*

**Automatic**  
measurements

**Panoramic view**  
of centimeter-level regions



## Why Image Analysis Method?

### Easy

Capture an image of particles, identify particles, then measure their size and shape. Every step of image analysis is easy and clear.

### Seeing is believing

The image analysis method determines the size and shape of every individual particle and then sums it up to form a statistic. Details of particle size or shape distribution can be accurately provided.

### Shape analysis

Based on a direct view of particles, it is possible to analyze not only the size of particles, but also their shape.

## Why Static Image Analysis Method?

### Clear vision

In static image analyzers, precision microscopes and cameras are specialized for high-quality particle images.

### Undersized particle sensitivity

The static image analysis method is sensitive to undersized particles; it is even possible to estimate the size of undersized particles.

### Small sample volume

The static image analysis method requires a small volume of samples. A few drops of emulsions or a few micrograms of powders are enough to do a measurement.

## BeVision Series: Precision in Particle Vision



### BeVision S1

Classical and versatile static image analyzer for wet and dry measurements.



### BeVision M1

Automated static image analyzer for wet and dry measurements.



### BeVision D2

Dynamic image analyzer for dry measurement.

	Static Image Analysis		Dynamic Image Analysis
	BeVision S1	BeVision M1	BeVision D2
Measurement range	1 - 3,000 μm	1 - 10,000 μm	30 - 10,000 μm
Particle shape analysis	●●●	●●●	●●●
High-resolution for narrow distributions	●●●	●●●	●●●
Accuracy for broad distributions	●	●●	●●●
Reproducibility	●	●●	●●●
Small sample volume for a single analysis	●●●	●●	●
Undersized particles detection	●●●	●●	●
Oversized particles detection	●	●	●●●
Simple operation and measurement efficiency	●●	●●●	●●●
Individual particle analysis	●●●	●●●	●●

# Efficient Scanning Mode and Limit - breaking Panoramic Mode



## Scanning Mode

The workflow of the BeVision M1 scanning mode is to capture an image first, then analyze the image while moving the stage, capture the next image once the stage has reached a new position, and repeat.

The BeVision software will display real-time results during the scanning process. The scanning mode is widely welcomed in different industries with its efficiency and reliability.

### Efficient and reliable scanning mode

Compared with the manual test, the automatic scanning process improves the test efficiency, doing the image capturing and stage moving simultaneously pushes the efficiency to the next level. The efficient scanning mode analyzes many particles in one test, thus strengthening the statistical significance of the result.




## Panoramic Mode

The panoramic mode is to stitch separate images into a panorama before analysis.

With a panorama, it is easy to measure the total number of particles in the region and the size and shape of every particle, even if the particle is an oversized one. The accurate particle number is of vital importance to different industries, e.g., the automotive industry and the pharmaceutical industry.

### A combination of macro vision with micro details

Panoramas produced by the BeVision M1 records all particles dispersed in a millimeter-level region and keeps their shape details. A clear panorama helps to accurately count particles in the macro-region and analyze the size and shape of particles, without any leftovers.



# Particle Size and Shape Parameters


## Size parameters

**Equivalent diameters:**  
area-equivalent diameter  
perimeter-equivalent diameter

**Feret diameters:**  
maximum and minimum Feret diameters,  $x_{LF}$  ("length")

**Martin diameters:**  
maximum and minimum Martin diameters

**Legendre ellipse:**  
major and minor axes




## Shape parameters

**Size difference in 2 directions:**  
aspect ratio  
L/W ratio  
ellipse ratio

**Round-likeness and rectangle-likeness:**  
circularity  
irregularity  
compactness  
extent  
box ratio

**Contour concavity:**  
concavity  
convexity  
solidity


**For elongated particles:**  
elongation  
straightness



# BeVision Software: Visualized Insights for You


## Reproducible measurements

To ensure a reproducible result, the BeVision software can make a measurement automatically, following a saved standard operation procedure (SOP).




## Comparable results

With the help of the BeVision software, it is possible to do a comparison among multiple records: particle size or shape distribution comparison, typical value comparison, etc.



## Customizable reports

The BeVision series prepares various report templates for different evaluation options. Layouts and contents of report templates are editable and customizable.




## Particle details

For irregularly shaped particles, it is hard to describe their size with a single dimension. Scanning over 180 different directions of each particle projection, the BeVision software is able to precisely analyze particles, and present the particle size and shape in 24 different parameters. The size and shape parameters are in compliance with ISO 9276-6.

## Distribution in total

Distribution curves and charts present particle size and shape distributions, and the scattering mode shows the relationship between two different particle size and shape parameters. All these charts, curves, and tables are customizable.



## Locate particles

The BeVision software offers a single particle gallery that can be the direct way to locate particles with a specific appearance. Besides, the BeVision software allows users to find particles with specific characteristics, with a customizable filter.


## Distribution in summary

The BeVision software offers statistics and typical values to describe particle size and shape distributions, e.g., the D[1,0], span value, and D90. The typical value chart is customizable.

# Application Cases


## Coolant tube cleanliness measurement

Surface cleanliness is a key factor for machinery manufacturing, the automotive industry, etc. As clarified in the ISO-16232 standard, an important step of coolant tube cleanliness measurement is to measure the size distribution of particles trapped on membrane filters. Also, the types of pollutants need to be identified, such as metals or microfibers. The BeVision M1 can do necessary inspection and analysis automatically, such as counting the particle number, identifying microfibers, and the cleanliness classification.




## Mineral pigments

The mineral pigment is a type of pigment derived from ground ore. The size and shape of mineral pigment particles affect the quality and performance of paint products. The BeVision M1 offers high-resolution size and shape measurement results at good efficiency, helping the QC engineers achieve an insightful validation of mineral pigments. A scatter plot showing the relationship between particle size and circularity helps compare the shape distribution of samples A and B, and guides designing paint recipes.



## Microorganisms

Traditionally, counting colonies on a petri dish requires taking a microscope image first, and counting the number then. A BeVision M1 could easily handle this job with the panoramic mode. Taking the panorama automatically, then measuring the total number of cells and analyzing the size and shape of them, everything will be done efficiently.



# Typical Applications

Agriculture



Abrasives



Paints, Inks & Coatings



Metal Powders



Surface Cleanliness



Mining and Minerals



Automotive



Ceramics



# BT - 910 Helps to Prepare Dry Powders



## How does it help?

The BT-910 dry powder dispersion module generates a pre-set air pressure difference, which drives the dispersion airflow. The BT-910 aims to offer a reliable and reproducible dispersion method for dry powders.

## Features and Benefits

- Reproducible dispersion
- No aggregates
- Even Dispersion

## General

Measuring principle	Static image analysis method, automatically scanning
Parameters	Particle size, shape, and number

## Measurement performance

Measuring range	1 – 10,000 $\mu\text{m}$
Typical measurement time	3 to 5 min *
Number of size/shape classes	100 (user adjustable)
Special functions	SOP settings, analysis of saved images

## Main device

Optical lens	4 x, 5 x, 10 x, 20 x (with 40 x digital magnification)
High-speed camera	Up to 120FPS
Light source	Halogen lamp, Köhler illumination
Stage scanning range	55 x 55 mm

## System parameters

Dimensions (L x W x H)	35.0 x 65.0 x 67.0 cm
Weight	18.7 kg
Supply voltage	100 / 240 V, 50 / 60 Hz

## Software

Conformity	ISO 13321, ISO 9276, ISO 16232, ISO 4406
Reports	Customizable reporting

\* Sample and sample preparation dependent

## BT – 910 dry powder dispersion module

Dimensions (L x W x H)	23.5 x 16.5 x 26.6 cm
Weight	4.3 kg
Supply voltage	100 / 240 V, 50 / 60 Hz
Dispersion air pressure	$\leq$ - 60 kPa

# Bettersize

BETTER PARTICLE SIZE SOLUTIONS

### Bettersize Instruments Ltd.

Website: <https://www.bettersizeinstruments.com>

Email: [info@bettersize.com](mailto:info@bettersize.com)

Address: No. 9, Ganquan Road, Jinquan Industrial Park,

Dandong, Liaoning, China

Postcode: 118009

Tel: +86-415-6163800

Fax: +86-415-6170645

Visit Our BeVision M1 Site:



Visit Our Official YouTube Channel:



**Disclaimer:** By using or accessing the brochure, you agree with the Disclaimer without any qualification or limitation. Diligent care has been used to ensure that the information in this brochure is accurate, Bettersize Instruments Ltd. shall not be liable for errors contained herein or for damages in connection with the use of this material. The information on this brochure is presented as general information and no representation or warranty is expressly or impliedly given as to its accuracy, completeness or correctness. It does not constitute part of a legal offer or contract. Bettersize Instruments Ltd. reserves the right to modify, alter, add and delete the content outlined in the brochure without prior notice and without any subsequent liability to the company.

Copyright: © 2022 Bettersize Instruments Ltd. | All Rights Reserved